





Portable Leaf Area Meter

By ADC BioScientific Ltd.



AM350 Portable Leaf Area Meter





- ✓ **SAVE TIME**. Do away with manual calculations for leaf area!
- ✓ You can rely on AM350 for accurate data, wherever your research takes you.
- ✓ Simply and easily view each live and stored image and data set.
- ✓ Name your files. Choose your averaging options for robust data.
- ✓ All data stored on internal memory.



What does AM350 Measure?





Area

Accumulated area

Mean area

Diseased leaf area

Length

Width

Perimeter

Ratio $R = L / W_m$

Shape factor $F = 4\pi A/P^2$

Units: Choice of mm, cm or inches







Anywhere

- ✓ Lightweight (1.8kg)
- ✓ Rechargeable battery (3,000 scans)
- ✓ Carrying bag
- ✓ Multi-position handle/stand

Anytime

- ✓ Large, real time image display
- ✓ High resolution (0.065mm²)
- ✓ High precision (1%)
- ✓ 2 year warranty

With Ease

- ✓ Menu driven software
- ✓ Vertical or horizontal mode
- ✓ Simple technique
- ✓ Image and data storage
- ✓ USB or serial download



AM350 Versatile Applications

- ✓ Whole leaf area on a variety of plants
- ✓ Diseased leaf and discoloured leaf area
- ✓ Insect damage
- ✓ Certain roots and rhizomorphs
- ✓ Extend scan using additional surface

for:

- ✓ Long leaves
- ✓ Wide leaves








Simple to use: Large Image Display





- ✓ View leaf scan in real time
- ✓ All parameters + image displayed
- ✓ High resolution image saved
- ✓ Choice of .bmp or .jpg files





Non-destructive Scans





- ✓ Clear, gentle membrane for positioning and smoothing leaf surfaces
- Hand scanner with smooth movement
- ✓ You have full control of amount of pressure applied to leaf
- ✓ Hold the lightweight AM350 up to leaves with shorter petioles, or, at canopy level





Measure Long and/or Wide Leaves

Adaptable to long and/or wide leaves:

For long, intact leaves up to 2m in length, the image screen will scroll to display and save entire leaf area, length and width (without perimeter and shape factor).

For wide leaves, scanning in 'strips' along the length will enable you to save each strip and add each to a running total. Total leaf area is automatically calculated.

Application Notes available from your local ADC support centre:

#AN001: Measuring Large, Wide Leaves with AM350 #AN010: Measurement of Diseased Areas of Leaves with AM350





#AN070: Measuring Long Leaves with AM350

www.hoskin.ca | sales@hoskin.ca Vancouver | Oakville | Montréal

Measure Diseased or Damaged Area





- ✓ Adjustable contrast control
- ✓ Measure:

Damaged area

Diseased area

Discoloured area

Application Note available:

"Measurement of diseased leaf area"



Portable: Adjustable Handle/Stand





Fully adjustable, locking carry handle and support stand





Makes AM350 easy to carry and position in varied locations



www.hoskin.ca | sales@hoskin.ca Vancouver | Oakville | Montréal

AM350 Protects your Data

Years of

Service to

Science BioScientific Ltd.

- ✓ Measurement in mm, cm, inches
- ✓ Data / image storage: 3000 measurements
- ✓ Data file management
- ✓ Data / image review on screen
- ✓ USB data download
- ✓ Data / image download in .bmp or .tif
- ✓ Automatic time and date storage
- ✓ Internal battery back up







AM350 Technical Specifications

Measured parameters: Leaf area, maximum length, maximum width,

perimeter, mean area, accumulated area, ratio and shape factor

Units of measurement: Select mm, cm, or inches

Scanner: Hand held unit with internal illuminator and CIS array.

Scanning speed: Up to 20mm/sec **Max. measurement width:** 103mm

Max. measuring length: 2m Precision / repeatability:

+/-1% Linear,

+/- 2% Area,

+/- 5% Perimeter

Resolution: 0.065mm2

Memory: 256KB RAM (around 2,000 data sets)

+ 560 scan images maximum (depending on complexity)



Display: 64 x 240 pixel LCD graphics

Battery: Nickel metal hydride 14.4V 1.6Ah

Internal battery back up

Computer interface: Mini-B USB connection and RS232

Battery charger: Built in fast charger; can use supplied mains adapter or

12V car battery.

Indicators for charge status. Charge time typically <3hrs.

Scans between charges: Typically 500+ scans from full charge

Operating temperature range: 0oC - 50oC

Dimensions of main unit: 275mm x 250mm x 30mm

Cable length: 90cm

Weight: 1.8kg (excluding charger)



WWW.HOSKIN.CA

• ENVIRONMENTAL • INSTRUMENTATION • MATERIALS TESTING

INTEGRATED SYSTEMS • RENTALS • SERVICE